



Hepatitis Profiles/Panels Testing Performed Monday - Friday	Antrim	Specimen Type/# of mL
HEPATITIS A PROFILE Tests: HAV Antibody (Total), HAV Antibody (IgM)	30060	Serum /1.0
HEPATITIS B PROFILE Tests: HBV Surface Antigen, HBV Surf Ab (Qual), HBV Core Ab (Total) Reactive HBV Core Ab (Total) will reflex HBC Core Ab (IgM)	31055	Serum /2.0
ACUTE HEPATITIS PANEL (ABC) Tests: HBV Surface Antigen, HBV Core Ab (IgM), HCV Antibody, HAV Antibody (IgM), HCV serum to MCM Reactive/Indeterminate HCV Antibody will reflex HCV RNA QN RT-PCR	32002	Serum/4.0
CHRONIC HEPATITIS PANEL (BC) Tests: HBV Surface Antigen, HBV Surf Ab (Qual), HBV Core Ab (Total), HCV Antibody, HCV serum to MCM Reactive/Indeterminate HCV Antibody will reflex HCV RNA QN RT-PCR	31065	Serum/4.0
ACUTE/CHRONIC HEPATITIS PANEL (ABC) Tests: HBV Surface Antigen, HBV Surf Ab (Qual), HBV Core Ab (Total), HBV Core Ab (IgM), HAV Antibody (IgM), HCV Antibody, HCV Serum toMCM Reactive/Indeterminate HCV Antibody will reflex HCV RNA QN RT-PCR	31060	Serum/5.0
NEEDLESTICK INCIDENT PROTOCOL – Source/Patient Tests: Rapid HIV-1/HIV-2 Antibody, HBV Surface Antigen, HCV Antibody Reactive/Indeterminate HCV Antibody will reflex HCV RNA QN RT- PCRReactive HIV-1/HIV-2 Antibody will reflex HIV Western Blot	30008	Serum/5.0 (Hepatitis) EDTA (HIV)
NEEDLESTICK INCIDENT PROTOCOL – Exposed/Employee Tests: Rapid HIV-1/HIV-2 Antibody, HBV Surf. Ab (Qual), HCV Antibody, HCV serum to MCM Reactive/Indeterminate HCV Antibody will reflex HCV RNA QN RT-PCR- Reactive HIV-1/HIV-2 Antibody will reflex HIV Western Blot	30009	Serum/5.0 (Hepatitis) EDTA (HIV)



METROPOLITAN
MEDICAL
LABORATORY, PLC

Individual Tests Testing Performed Monday - Friday	Antrim	Specimen Type/# of mL
HAV Antibody (Total) Alias: anti-HAV (Total), Hepatitis A IgG & IgM Will reflex to anti-HAV-IgM if positive.	30025	Serum /1.0
HAV Antibody (IgM) Alias: anti-HAV (IgM)	30030	Serum /1.0
HBV Surface Antigen Alias: HBsAg Reactive HBV Surface Antigen will reflex HBV Surf Ag Confirmatory Assay	30000	Serum /1.5
HBV Surface Antibody (Qualitative) Alias: anti-HBS	30015	Serum /1.0
HBV Surface Antibody (Quantitative) Alias: anti-HBS (Quant), QHBSAB	30040	Serum /1.0
HBV Core Antibody (Total) Alias: ANTI-HBC, ANTI-HBV CORE Will reflex to Anti-HBc-IgM if positive.	30020	Serum/0.5
HBV Core Antibody (IgM) Alias: Anti-HBC IgM	30035	Serum /1.0
HBV Core Antibody (Total) Alias: anti-HBc Total (no reflex)	30021	Serum/1.0
HCV Antibody , HCV Serum to MCM Alias: anti-HCV Reactive/Indeterminate HCV Antibody will w/reflex to HCV RNA QN RT-PCR	32080	Serum /3.0
HIV-1/HIV-2 Antibody Alias: anti-HIV-1/-2 Reactive HIV-1/HIV-2 Antibody will reflex HIV Western Blot	20020	Serum /1.0
HIV-1/HIV-2 Antibody (Rapid) Alias: anti-HIV-1/-2 Reactive HIV-1/HIV-2 Antibody will reflex HIV Western Blot	106350	EDTA



METROPOLITAN
MEDICAL
LABORATORY, PLC

Referred Tests	Antrim	Specimen Type/ # of mL
Hepatitis Be Antigen (HBeAg)	302135	Serum /1.0
Hepatitis Be Antibody (anti-Hbe)	302130	Serum /1.0
HCV RNA Detection/Quantitation ONLY	306061	Serum/Plasma3 .0
Hepatitis Delta Virus Antibody (anti-HDV)	302345	Serum or Plasma /2.0
Hepatitis Delta Virus (HDV), Antigen	303260	Serum/2.0
HIV-1/HIV-2 Confirm Eval (aka: HIV Western Blot)	304734	Serum/1.0
HBV DNA Quantitative by RT-PCR (aka: HBV Monitor, HBV viral load)	303120	Serum or Plasma/3.0
HIV-1 RNA by RT-PCR Quantitative (aka: HIV monitor, HIV viral load)	302565	EDTA PLASMA/ 2.0
HBV Genotyping	306005	Serum or Plasma/1.0
HCV Genotyping	302465	SerumorPlasma /2.0
HIV-1 Genotyping	304365	EDTA PLASMA/ 2.0